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1) H. Matsushita, Y. Takahashi, S. Azuma, H. Hiroi and T. Amagai: Development of highly sensitive automatic analysis for polynuclear aromatic hydrocarbons in airborne particulates and its application to the survey of indoor pollution. Indoor Air Quality in Asia (G. B. Leslie, K. J. Leslie, J. Huang, Y. Qin), pp236-243, Indoor Air International, International Association for Indoor Air Quality, Rothenfluh Switzerland (1994).

2) Lloyd R. Snyder, Joseph L. Glajch and Joseph J. Kirkland: Practical HPLC Method Development. John Wiley & Sons, Inc, 1988. (訳本: 高速液体クロマトグラフィーの実際、高橋昭、荒木峻訳、東京化学同人。)


3) US EPA Method To-13, "The determination of Benzo[a]pyrene (BaP) and other polynuclear aromatic hydrocarbons (PAH's) in ambient air using gas chromatographic (GC) and high performance liquid chromatographic (HPLC) analysis". 1988.

Procedure for measurement of ambient NO₂ concentration by the Saltzman's method


Principle

This method is originally developed by B.E. Saltzman (1954, Anal. Chem., 36, 610) and modified several times, especially about conversion ratio of NO₂ gas to nitrite ion. Now conversion ratio (Saltzman Factor) is experimentally determined as 0.84. Namely, 1 mol NO₂ produces 0.84 mol NO₂⁻ ion in a Saltzman solution. The absorbing solution includes sulfanilic acid (p-amino benzenesulfonic acid) and N-(1-naphtyl) ethylenediamine in 5 % acetic acid solution.

The NO₂ in a sample air collected in a Saltzman solution reacts with sulfanilic acid and produces diazonium compound (diazotation reaction).

Followed the above reaction, the compound reacts with N-(1-naphtyl)ethylenediamine (NEDA) and produces an azo dye which is stable red color (coupling reaction). This colored solution is measured by spectrophotometry.

Reagent

1. Absorbing solution (durable 3 months)

Sulfanilic acid (SA)	5.0 g	
N-(1-Naphtyl) ethylenediamine dihydrochloride (NEDA)		0.05 g
Acetic Acid	50 ml	
Distilled water	950 ml	

SA and NEDA are dissolved completely in 950 ml distilled water using a glass stick, then add 50 ml acetic acid. This solution is transferred to a 1 l brown bottle and stored in a fridge.

2. NO₂ standard solution (durable 1 day)

Dry NaNO₂ at 110 °C for three hours, then keep in a silicagel desiccator for one night.

Weigh 0.25 g and dissolved in 1 l of distilled water (original solution).

Take 1 ml of original solution in a 100 ml volumetric flask and fill distilled water until the line. This final solution (standard solution) 1 ml equivalents to 1.0 µl NO₂ gas (0 °C, 760 mmHg)

NO₂ concentration in original solution is calculated as follows.

$$\frac{0.25 \text{ g}}{69 \text{ g}} \times 22.4 \text{ l} \times \frac{1}{1000 \text{ ml}} \times \frac{1}{0.84} = 1 \times 10^{-4} \text{ l/ml}$$

The standard solution must be prepared just before use, however, the original solution can be stored in a refrigerator for one week.

Analytical Line

Take 1, 2, 3, 4, 5 ml of the standard solution in 20 ml of graduated test tubes each, then, add 9, 8, 7, 6, 5 ml of absorbing solution and leave 20 minutes at room temperature. Then transfer to photo-cell tubes for the spectrophotometer and measure the absorbance at wave length 540 nm.

Each absorbance is plotted on a normal scale section paper with the vertical axis for absorbance and the horizontal for amount of NO₂.

In case the points are deviate from the regression line, this work should be repeated.

Procedure for measuring absorbance

A simple type of spectrophotometer is shown in Fig. 1.


Fig.1. A simple type of Spectrophotometer

Turn the dial B clockwise to power on state

Turn the dial D and adjust wave length at 540 nm.

Turn the dial A and adjust indicator to ∞ then open lid C and insert a blank tube with absorbing solution and adjust the indicator to 0 using dial B. Repeat above procedure several times and confirm the stability. Measure the absorbance of sample solutions. (Notice: Absorbance are written in red figures).

Procedure of sampling

Ten ml of absorbing solution is taken in an absorbing tube (brown color) and connect with an empty impinger and a handy air sampler as shown in Fig.2.

Push the power switch A and adjust the flow rate 0.3 l/min by the regulator B.


Fig.2. Schematic diagram of the apparatus for measuring NO₂ by Saltzman's method

Calculation

$$C_{\text{NO}_2} (\text{ppm}) = \frac{A}{A_s} \times \frac{1}{v}$$

C_{NO_2} (ppm) : NO₂ concentration (ppm)

A : Absorbance of a sample in 10 ml solution

A_s : Absorbance of standard solution contains 1 μl (NTP) of NO₂ in 10 ml of solution (nearly 0.19)

v : Volume of sampled gas (1, 0 °C, 1 atm)