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Fast Analysis of Hydrophilic Compounds in Environmental Waters by use of activated carbon fiber cartridge.

Masaru Yasuda, Kuniaki Shimbo, Akiko Shinoda and Hiroshi Suzuki
Showa Denko K.K.

Solid phase extraction cartridges (SPE), such as C18 or activated carbon, have been utilized for the sample preparation for many compounds in environmental waters; nowadays demands of many researchers are to make up faster sample preparation.

Thus, we have been trying to develop a new SPE cartridge packed with activated carbon fiber (AC-F) to achieve faster analysis for many compounds in environmental waters compared with the conventional activated carbon cartridge.

The AC-F cartridge can be used at 300mL/min as a flow rate after passing through 800mL of 0.001% gelatin containing water. After 1L river water containing 20 μ g of methomyl is passed through the AC-F cartridge, methomyl is eluted from the cartridge by organic solvent and checked its recovery by HPLC. When the flow rates are 32mL/min and 50mL/min, the recoveries of methomyl are 80-90% and 80% or more, respectively.

In this study, we tried to check the recovery of 1,4-dioxane, which is a typical hydrophilic compound in purified water. When flow rate are 20mL/min and 40mL/min, the recoveries of 1,4-dioxane are 84% and 80% or more, respectively.