
preserved buffer.

During the assay procedure, magnetic separation was performed using a magnetic separation rack (Abraxis, Warminster, PA). The device consists of a magnetic separation rack and a test tube holder, which fits over a magnetic separation rack containing permanent rare earth magnets. The two-piece design allows for up to sixty tubes to be set up, incubated and magnetically separated without removing the tubes from the tube holder.

Water samples were analyzed directly without any sample extraction or pre-concentration. After mixing 1:1 with methanol, 250 mL of the sample and anti-PBDE coupled magnetic particles (500 mL) were added to a disposable glass test tube and incubated for 20 minutes at room temperature. An aliquot of PBDE-HRP solution (250 mL) was added and the reaction was incubated for 20 minutes at room temperature. A magnetic field was applied to the magnetic solid-phase to facilitate washing and removal of unbound PBDE-HRP and to eliminate any potential interfering substances. The enzyme substrate and chromogen (peroxide/TMB) were then added (500 mL) and incubated for 20 minutes. The reaction is stopped with the addition of 2 N H<sub>2</sub>SO<sub>4</sub> (500 mL) and the final color was read with analyzed using a photometer by determining the absorbance at 450 nm. The observed absorbance data were compared to a linear regression line using a log-log standard curve prepared from calibrators containing 0, 25, 50, 100, 500 and 1,000 parts per trillion (ppt) of BDE-47.

**Results and Discussion**

*Dose Response Curve and Sensitivity.* Figure 1 illustrates the mean standard curve for the BDE-47 calibrators collected over 24 assays, error bars represent one standard deviation (SD, n=24). The displacement at the 25 ppt level is significant (84% B/Bo, where B/Bo is the absorbance at 450 nm observed for a sample or standard divided by the absorbance at the zero standard). The assay sensitivity in water based on 90% B/Bo is 17 ppt.

*Precision.* Table 1 shows the results from a precision study in which surface and groundwater samples fortified with BDE-47 at 3 concentrations were each measured 5 times per assay and on five different days. The within and between day variation was estimated by the method of Bookbinder and Panosian<sup>13</sup>. Coefficients of variation below 5% were observed.

**Table 1 Table 2**

**Precision of BDE-47 Measurements in Water BDE-47 Recovery in Water**

Conc. of BDE-47 Added (ppt)	Recovery		
	Mean (ppt)	S.D. (ppt)	%
62.5	56.4	2.9	90
125	131.5	5.3	105
250	270.6	12.8	108
500	510.4	28.3	102
Average			101.4

Control	1	2	3
Replicates	5	5	5
Days	5	5	5
n	25	25	25
Mean (ppt)	47	95	514
% CV (within assay)	3.2	2.2	2.0
% CV (between assay)	4.4	4.8	4.0

*Accuracy.* Known amounts of BDE-47 were added to four groundwater samples obtained from Warminster, PA. The samples included a municipal water source, a reservoir, a lake, a pond,

and a creek. The accuracy (recovery) was assessed by analyzing the samples before and after the addition of BDE-47 and then subtracting the estimated concentration of PBDE before spiking. None of the samples had significant levels of PBDEs. Added amounts were accurately recovered (Table 2). An average assay recovery of 101% was obtained.

**Table 3**

**Specificity (Cross-Reactivity)**

Compound	LDD (ppb)	50% B/Bo (ppb)
BDE -47	0.017	0.135
BDE-99	0.02	0.15
BDE -28	0.045	0.9
BDE 100	0.055	5.5
BDE 153	0.075	10
BDE 154	3.5	580
BDE -183	13.5	2000
BDE -209	370	3000
Aroclor 1254	3	180
PCB-37	160	2000
PCB-77	>1000	>1000
PCP	3300	>10,000
2,4-D	>10,000	>10,000

**Specificity.** Table 3 summarizes the cross-reactivity data of the PBDE assay for various PBDE congeners as well as other environmental contaminants such as PCBs, PCP, and 2,4-D. The percent cross-reactivity was determined by estimating the amount of analogue required to displace 50% of the enzyme conjugate and comparing to the 50% displacement of the BDE-47 standard curve. The results showed that the antibody recognizes BDE-47 and BDE-99 equally well, the higher the brominate substitution the less cross-reactivity was observed. The antibodies have less than 0.1% cross-reactivity for Aroclor 1254.

### Conclusions

This work describes the performance characteristics of a magnetic particle-based ELISA for the detection of PBDEs in groundwater samples. The assay is fast, and eliminates the need for expensive instrumentation and solvent disposal. The ELISA exhibits good precision and accuracy which can provide consistent and cost-effective monitoring of water

samples. Using this ELISA, fifty results can be obtained in about one hour. Future efforts will be to extend these observations to the analysis of other type of matrices such as food, serum, milk, soil, and to perform comparison with other analytical methods.

### Acknowledgements

The authors would like to thank Amy McGarvey for technical assistant.

### References

1. IPCS (1994). Environmental Health Criteria 162: Brominated diphenyl ethers. International Programme on Chemical Safety, World Health Organization, Geneva.
2. Hardy, M.I. *Chemosphere* **2002**, *46*, 757-777.
3. deWitt, C.A. *Chemosphere* **2002**, *46*, 583-624.
4. Boon, J.P.; Lewis, W.E.; Tjoen-A-Choy, M.R.; Allchin, C.B.; Law, R.J.; De Boer, J.; Hallers-Tjabbes, C.C.; Zegers, B.N.; *Environ. Sci. Technol.* **2002**, *46*, 697-707.
5. Anderson, O., and Blomkvist, G. *Chemosphere* **1981**, *10*, 1051-1060.
6. Klasson Wehler, E.; Hovander, L.; Bergman, A. *Organohalogen Compd.* **1997**, *33*, 420-425.
7. Meiroyte, D.; Noren, K; Bergman, A. *J. Toxicol. Environ. Health.* **1999**, *58*, 329-341.
8. Brouwer, A.; Morse, D.C.; Lnas, M.C.; Schuur, A.G.; Murk, A.J.; Klasson-Wehler, E.; Bergman, A., and Visser, T.J. **1998**, *Toxicol. Ind. Health*, *14*, 59-84.
9. Hottenstein, C.S.; Fleeker, J.R; Herzog, D.P.; Rubio, F.M.; Lawruk, T.S. *Journal of Agricultural and Food Chemistry* **1996**, *44*:3576-358.
10. Jourdan, S.W.; Scutellaro, A.M; Fleeker, J.R.; Herzog, D.P.; Rubio, F.M. *Journal of Agricultural and Food Chemistry* **1995**, *43*:2784-2788.
11. Lawruk, T.S.; Gueco, A.M.; Jourdan, S.W.; Scutellaro, A.M.; Fleeker, J.R.; Herzog, D.P.; Rubio, F.M. *Journal of Wine Research* **1994**, *5*:205-214.
12. Shelver, W. L., Keum, Y. -S., Kim, H. -J., Rutherford, D., Hakk, H., Bergman, D., and Li, Q. X. *J. Agri. Food*

ANA - Bioanalytical Approaches for POPs Detection

*Chem.* **2005** (web release: April 15, 2005)

13. Bookbinder, M.J.; Panosian, K.J. Correct and Incorrect Estimation of Within and Between-day Variation. *Clin. Chem.* **1986**, *32*:1734-1737.

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#### 抗PBDE抗体固相化磁気微粒子(PBDE Antibody Coupled Paramagnetic Particles)

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